

**Centrifugal stirring rotors, 2-blade, stainless steel 1.4305****Centrifugal stirring rotors, 2-blade, stainless steel 1.4305**

The folding blades of this 2-blade centrifugal stirring rotor open with increasing speed, allowing use with narrow-neck vessels. The stirring effect it provides is similar to a standard diagonal blade stirring rotor. Recommended for medium to high speed operations.



Code	Description	Packaging
LLG09197011	Centrifugal stirring rotors, 2-blade, stainless steel 1.4305, Ø agitator open / closed 60/15 mm, Ø agitator shaft 8 mm, Length 350 mm	1 pz.
LLG09197013	Centrifugal stirring rotors, 2-blade, stainless steel 1.4305, Ø agitator open / closed 100/24 mm, Ø agitator shaft 8 mm, Length 550 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

