



Flask support rings, "BiBase", silicone elastomer

Flask support rings, "BiBase", silicone elastomer

Blue. Double sided. Sits firmly on the bench to support round-bottomed flasks.

- Resistant to oxidation and temperatures from -60 to 200°C
- Very soft and easy to clean
- Resistant to chemicals
- Good thermal and electric insulation
- Excellent adhesion with glass
- Non-slip



Code	Description	Packaging
LLG09143035	Flask support rings, "BiBase", silicone elastomer, For flasks 50 ... 500 ml, Ext. diam. 80 mm, Int. diam. 46 mm, Height 30 mm, Max. 13 kg	1 pz.
LLG09143037	Flask support rings, "BiBase", silicone elastomer, For flasks 2000 ... 6000 ml, Ext. diam. 162 mm, Int. diam. 120 mm, Height 45 mm, Max. 30 kg	1 pz.
LLG09143036	Flask support rings, "BiBase", silicone elastomer, For flasks 500 ... 2000 ml, Ext. diam. 120 mm, Int. diam. 80 mm, Height 40 mm, Max. 18 kg	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

