

**Universal centrifuge Z 366 / Z 366 K****Hermle****Universal centrifuge Z 366 / Z 366 K**

These universal table top centrifuges Z 366 and Z 366 K putting new scales regarding the versatility, running features, easy handling and design. Both, the tube capacity of 40 x 15 ml conical tubes and centrifugal forces up to 24325 xg, turning these table top centrifuges Z 366 and Z 366 K into universal high-powered centrifuges.

- Maintenance-free induction drive
- Microprocessor with large LCD Display
- Motor driven lid lock
- Automatic rotor identification system with over speed protection, immediately upon rotor insertion
- Active imbalance-identification and cut off
- CFC-free refrigeration system
- A large variety of rotors and accessories
- Manufactured according to international safety regulations, i.e. IEC 61010
- Audible signal at the end of each run with several melody options
- Rotor exchange is easy to handle
- Noise level 63 dBA (depends on rotor)
- Optimized airflow system keeps the samples cool

Without rotor. Please order separately.

Specifications

Max. RCF Z 366 / Z 366 K: 21379 / 24325 x g

Max. speed Z 366 / Z 366 K: 15000 / 16000 rpm

Max. capacity: 6 x 250 ml

Cooling temperature Z 366 K: -20 °C to 40 °C

Timer: 59 min 59 sec / 10 sec increments

99 h 59 min / 1 min increments

Power consumption Z 366 / Z 366 K: 530 W / 800 W

Dimensions (W x D x H) Z 366 / Z 366 K: 430 x 510 x 360 mm / 720 x 510 x 360 mm

Weight Z 366 / Z 366 K: 52 kg / 77 kg

Power supply: 230 V, 50/60 Hz



Code	Description	Packaging
LLG06268605	Universal centrifuge Z 366 / Z 366 K, Type Universal Centrifuge Z 366 K	1 pz.
LLG06269721	Universal centrifuge Z 366 / Z 366 K, Type Universal Centrifuge Z 366, 120 V, 50/60 Hz	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

