

**Parts for Spectrophotometers, 6700 series****Jenway****Parts for Spectrophotometers, 6700 series****A range of three Split-Beam instruments:**

Type 6700: Visible range scanning (320nm to 1100nm with 4nm spectral bandwidth)

Type 6705: UV/Visible range scanning (190nm to 1100nm with 4nm spectral bandwidth)

Type 6715: UV/Visible range scanning (190nm to 1100nm with 1.5nm spectral bandwidth)

Modes of operation: Photometrics, Kinetics, Quantitation, Multiwavelength analysis and Spectral scanning.

Featuring:

- Split-Beam Spectrophotometer
- Fast scanning rate of 1500nm/minute with data collection at 0.1nm intervals
- Up to 20 point calibration for Quantitative analysis
- Multiwavelength mode for DNA/RNA analysis
- Secure Multi-User Operation, full password protection
- Colour (TFT) LCD with Touch-Screen Interface and sealed QWheeläControl
- Full post scan analysis (Peak / Valley pick, Peak Ratios, Area, Zoom, Wavelength table, Derivatives, Smoothing)
- Removable SD / Multi-Media Memory Card, up to 2GB available
- Plug-In Accessory Modules including 8 cell changer, water heated 6 cell changer, Peltier devices for temperature controlled measurements and sipper systems for high throughput analysis. A full range of cuvette and test tube holders are also available.
- 3 year warranty



Free PC software supplied for data transfer.

Supplied with: Mains lead, pack of 100 disposable cuvettes, 10 x 10mm cell holder, PC Application Software on CD-ROM and operating instructions.



Code	Description	Packaging
LLG09775434	Parts for Spectrophotometers, 6700 series, Type for all types, Printer, incl. 1 roll of paper	1 pz.
LLG06253506	Parts for Spectrophotometers, 6700 series, Type for all types, Sample chamber with automatic 8-cell	1 pz.
LLG07900098	Parts for Spectrophotometers, 6700 series, Type for all types, Sample chamber with Peltier module fitted	1 pz.
LLG09775435	Parts for Spectrophotometers, 6700 series, Type for all types, Sample chamber with sipper module fitted	1 pz.
LLG06236533	Parts for Spectrophotometers, 6700 series, Type for all types, Sample chamber with sipper Peltier module fitted	1 pz.
LLG06236532	Parts for Spectrophotometers, 6700 series, Type for all types, Sample chamber with water heated automatic six- cell changer module fitted	1 pz.
LLG09775436	Parts for Spectrophotometers, 6700 series, Type for all types, Test tube holder for 10mm square cuvettes, 16mm and 24mm diameter tubes	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

