


**LLG-Draining racks, PE-coated wire****LLG Labware****LLG-Draining racks, PE-coated wire**

Wall or bench mounting for 50 reagent bottles and 5 flasks

With drip channel and plastic covered feet.

Dimensions (L x W x H): 360 x 130 x 650 mm.



Code	Description	Packaging
LLG06281619	Draining racks, PE-coated wire, For 50 reagent bottles and 5 flasks	1 pz
LLG09041619	pH plastic electrode LE407with gel elec	1 pz.
LLG09041620	pH glass elektrode LE409with liquid ele	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

