

**Metal-cased heating mantles series KM-M, without controller****SAF Wärmetechnik****Metal-cased heating mantles series KM-M, without controller**

LabHEAT® - Heating mantles for flasks, round bottom. Flexible glass yarn heating element in plastic coated, chemical-resistant metal housing with built-in power-on and heating zone switch; max. heating element temperature 450°C; 1.5m power supply (earthed) cable and RCD (residual current detection); nominal voltage 230V AC.



Code	Description	Packaging
LLG09642509	Metal-cased heating mantles series KM-M, without controller, 6000 ml, Flask diam. 236 mm, Rating 1000 W, Heating zones 2	1 pz.
LLG09642501	Metal-cased heating mantles series KM-M, without controller, 100 ml, Flask diam. 64 mm, Rating 100 W, Heating zones 1	1 pz.
LLG09642504	Metal-cased heating mantles series KM-M, without controller, 1000 ml, Flask diam. 131 mm, Rating 300 W, Heating zones 2	1 pz.
LLG09642510	Metal-cased heating mantles series KM-M, without controller, 10000 ml, Flask diam. 279 mm, Rating 1400 W, Heating zones 2	1 pz.
LLG09642505	Metal-cased heating mantles series KM-M, without controller, 2000 ml, Flask diam. 166 mm, Rating 500 W, Heating zones 2	1 pz.
LLG09642511	Metal-cased heating mantles series KM-M, without controller, 20000 ml, Flask diam. 345 mm, Rating 2000 W, Heating zones 2	1 pz.
LLG09642502	Metal-cased heating mantles series KM-M, without controller, 250 ml, Flask diam. 85 mm, Rating 150 W, Heating zones 2	1 pz.
LLG09642506	Metal-cased heating mantles series KM-M, without controller, 3000 ml, Flask diam. 185 mm, Rating 600 W, Heating zones 2	1 pz.
LLG09642507	Metal-cased heating mantles series KM-M, without controller, 4000 ml, Flask diam. 207 mm, Rating 750 W, Heating zones 2	1 pz.
LLG09642500	Metal-cased heating mantles series KM-M, without controller, 50 ml, Flask diam. 51 mm, Rating 55 W, Heating zones 1	1 pz.
LLG09642503	Metal-cased heating mantles series KM-M, without controller, 500 ml, Flask diam. 105 mm, Rating 200 W, Heating zones 2	1 pz.
LLG09642508	Metal-cased heating mantles series KM-M, without controller, 5000 ml, Flask diam. 223 mm, Rating 860 W, Heating zones 2	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

