



Temperature sensors for Laboratory regulator series KM-RX

Temperature sensors for Laboratory regulator series KM-RX


LabHEAT®-sensors. Suitable for temperature control and KM-KM-RX1001
RX1004 with permanently attached cable, with or without diode plug.

KM-TP2: Pt100, up to 400 °C, sheated element V2A, dia. 4 x 50 mm,

KM-TPG: Pt100, up to 250 °C, glass tube, dia. 6 x 40 0mm,

KM-TNS: NiCr-Ni, up to 1200 °C, rod sensor in Inconell, dia. 1.5 x 300 mm

KM-TNF: NiCr-Ni, up to 400 °C, flat sensor, 0.4 x 15 x 400 mm



Code	Description	Packaging
LLG06236355	Temperature sensors for Laboratory regulator series KM-RX, Type KM-TNF, For KM-RX1001, Length 1.5 m, with diode plug	1 pz.
LLG06237998	Temperature sensors for Laboratory regulator series KM-RX, Type KM-TNS, For KM-RX1001, Length 2.0 m, with diode plug	1 pz.
LLG09642674	Temperature sensors for Laboratory regulator series KM-RX, Type KM-TPG, For KM-RX1001, Length 1.5 m, with diode plug	1 pz.
LLG09642670	Temperature sensors for Laboratory regulator series KM-RX, Type KM-TNF, For KM-RX1004, Length 1.5 m, without diode plug	1 pz.
LLG09642671	Temperature sensors for Laboratory regulator series KM-RX, Type KM-TNS, For KM-RX1004, Length 2.0 m, without diode plug	1 pz.
LLG09642672	Temperature sensors for Laboratory regulator series KM-RX, Type KM-TP2, For KM-RX1004, Length 2.0 m, without diode plug	1 pz.
LLG09642673	Temperature sensors for Laboratory regulator series KM-RX, Type KM-TPG, For KM-RX1004, Length 1.5 m, without diode plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

