


**Support clamp KM-SK****SAF Wärmetechnik****Support clamp KM-SK**

LabHEAT®-Accessories for metal-cased heating mantles and serial heating unit to fix support rods up to diameter 12 mm or to integrate the metal-cased heating mantles within support wall.



| Code | Description | Packaging |
|-------------|---------------------------------|-----------|
| LLG09642630 | Support clamp KM-SK, Type KM-SK | 1 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

