

**Glass fibre-insulated heating tapes series KM-HT-BS30****SAF Wärmetechnik****Glass fibre-insulated heating tapes series KM-HT-BS30**

LabHEAT®-Heating tapes. No protection against ingress of water; with inner metal protective braiding and glass fibre insulation; product dimensions 30 x 5mm; min. bending radius > 15mm; max. heating element temperature 450°C; preterminated with 1.0m cold end; nominal voltage 230V AC.



Code	Description	Packaging
LLG09642726	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 400 cm, Heating power 1000 W	1 pz.
LLG09642720	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 50 cm, Heating power 125 W	1 pz.
LLG09642729	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 1000 cm, Heating power 2000 W	1 pz.
LLG09642721	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 100 cm, Heating power 250 W	1 pz.
LLG09642722	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 150 cm, Heating power 375 W	1 pz.
LLG09642723	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 200 cm, Heating power 500 W	1 pz.
LLG09642724	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 250 cm, Heating power 625 W	1 pz.
LLG09642725	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 300 cm, Heating power 750 W	1 pz.
LLG09642727	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 500 cm, Heating power 1250 W	1 pz.
LLG09642728	Glass fibre-insulated heating tapes series KM-HT-BS30, Length 700 cm, Heating power 1500 W	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

