


**Glass fibre-insulated heating cables series KM-HC-G****SAF Wärmetechnik****Glass fibre-insulated heating cables series KM-HC-G**

LabHEAT®-Heating cables. No protection against ingress of water; glass fibre insulation without protective braiding; outer diameters dia. 3.5 to 4.5mm; min. bending radius > 5mm; max. heating element temperature 450°C; preterminated with 2 x 1.5m cold end; nominal voltage 230V AC.



Code	Description	Packaging
LLG09642751	Glass fibre-insulated heating cables series KM-HC-G, Length 100 cm, Heating power 150 W	1 pz.
LLG09642752	Glass fibre-insulated heating cables series KM-HC-G, Length 150 cm, Heating power 225 W	1 pz.
LLG09642753	Glass fibre-insulated heating cables series KM-HC-G, Length 200 cm, Heating power 300 W	1 pz.
LLG09642754	Glass fibre-insulated heating cables series KM-HC-G, Length 250 cm, Heating power 375 W	1 pz.
LLG09642755	Glass fibre-insulated heating cables series KM-HC-G, Length 300 cm, Heating power 450 W	1 pz.
LLG09642756	Glass fibre-insulated heating cables series KM-HC-G, Length 400 cm, Heating power 600 W	1 pz.
LLG09642750	Glass fibre-insulated heating cables series KM-HC-G, Length 50 cm, Heating power 75 W	1 pz.
LLG09642757	Glass fibre-insulated heating cables series KM-HC-G, Length 500 cm, Heating power 750 W	1 pz.
LLG09642758	Glass fibre-insulated heating cables series KM-HC-G, Length 600 cm, Heating power 900 W	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

