

**pH buffer solutions, technical**

Without colour coding, hot steam sterilized. Stated pH values at 25°C. Available in 250ml containers. Shelf life approx. 2 years, accuracy ± 0.02 pH.



Code	Description	Packaging
LLG09041043	pH buffer solutions, technical, pHvalue 4.00 at 25 °C, Colour Clear, 250 ml	1 pz.
LLG09041045	pH buffer solutions, technical, pHvalue 7.00 at 25 °C, Colour Clear, 250 ml	1 pz.
LLG09041048	pH buffer solutions, technical, pHvalue 10.01 at 25 °C, Colour Clear, 250 ml	1 pz.
LLG09041081	pH buffer solutions in ampoules, technical, pHvalue 4.00 / 7.00 / 10.01 at 25 °C, Package contents 3 x 20 ampoules	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

