



Microscope slide boxes

Microscope slide boxes

Colour-coded for easy identification. Slide boxes are perfect for storage or transportation of standard size slides (76mm x 26mm). With cork lining. Stackable. Each 100 place box is secured with a rust-resistant nickel-plated clasp and hinge pin. The 50 place box is secured by a thumb-latch lock mechanism. Each slot is numbered to correspond to the slide inventory sheet on the inside cover of the box.



Code	Description	Packaging
LLG09161287	Microscope slide boxes, To hold microscope 100 slides, Colour Black, dim. (W x D x H) 208 x 175 x 34 mm	1 pz.
LLG06241030	Microscope slide boxes, To hold microscope 25 slides, Colour Blue, dim. (W x D x H) 141 x 88 x 35 mm	1 pz.
LLG09161284	Microscope slide boxes, To hold microscope 100 slides, Colour Blue, dim. (W x D x H) 208 x 175 x 34 mm	1 pz.
LLG09161300	Microscope slide boxes, To hold microscope 50 slides, Colour Blue, dim. (W x D x H) 209 x 86 x 35 mm	1 pz.
LLG06240819	Microscope slide boxes, To hold microscope 25 slides, Colour Green, dim. (W x D x H) 141 x 88 x 35 mm	1 pz.
LLG09161285	Microscope slide boxes, To hold microscope 100 slides, Colour Green, dim. (W x D x H) 208 x 175 x 34 mm	1 pz.
LLG06240818	Microscope slide boxes, To hold microscope 25 slides, Colour Red, dim. (W x D x H) 141 x 88 x 35 mm	1 pz.
LLG09161286	Microscope slide boxes, To hold microscope 100 slides, Colour Red, dim. (W x D x H) 208 x 175 x 34 mm	1 pz.
LLG09161289	Microscope slide boxes, To hold microscope 50 slides, Colour Red, dim. (W x D x H) 209 x 86 x 35 mm	1 pz.
LLG06268500	Microscope slide boxes, To hold microscope 50 slides, Colour White, dim. (W x D x H) 209 x 86 x 35 mm	1 pz.
LLG07079563	Microscope slide boxes, To hold microscope 100 slides, Colour White, dim. (W x D x H) 208 x 175 x 34 mm	1 pz.
LLG09161288	Microscope slide boxes, To hold microscope 100 slides, Colour Yellow, dim. (W x D x H) 208 x 175 x 34 mm	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

