

**Wide-mouth Bottle Nalgene PassPort IP2, HDPE, with screw cap, PP****Thermo Scientific****Wide-mouth Bottle Nalgeneä PassPort IP2, HDPE, with screw cap, PP**

Safe, durable, cost-efficient line of wide-mouth bottles for packaging and transporting hazardous materials. Recommended for customers who are designing, assembling and certifying their own combination packaging. Split and puncture resistant. Wide mouth allows for easy filling of dry materials or liquids. Heavy duty walls. Moulded into base: IP2 marking, material, volume, registration notch and the Nalgene® brand name. Evaluated at 15psi (103kPa) per 49 CFR 173.27 (c) (2), ICAO Technical Instructions Part 3-1;1.1.6.1, and IATA Dangerous Goods Regulations Section 5.0.2.9. Leakproof.



Code	Description	Packaging
LLG09103247	Wide-mouth Bottle Nalgene PassPort IP2, HDPE, with screw cap, PP, 1000 ml	6 pz.
LLG09103244	Wide-mouth Bottle Nalgene PassPort IP2, HDPE, with screw cap, PP, 125 ml	12 pz.
LLG09103245	Wide-mouth Bottle Nalgene PassPort IP2, HDPE, with screw cap, PP, 250 ml	12 pz.
LLG09103246	Wide-mouth Bottle Nalgene PassPort IP2, HDPE, with screw cap, PP, 500 ml	12 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

