


**Stopcocks, with glass plug, borosilicate glass 3.3****Lenz****Stopcocks, with glass plug, borosilicate glass 3.3**

Acc to DIN 12541-1 and in addition to DIN standard. With solid glass plug and retaining device. Straight bore. Leakage rate acc. to DIN 12540-1.



Code	Description	Packaging
LLG07018947	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 12.5, Bore 1.5 mm, DIN, Material SBW	1 pz.
LLG09115611	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 12.5, Bore 1.5 mm, DIN 12541-1, Material	1 pz.
LLG06079648	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 14.5, Bore 2.5 mm, DIN 12541-1, Material	1 pz.
LLG07018936	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 12.5, Bore 2.5 mm, DIN, Material SBW	1 pz.
LLG07018948	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 14.5, Bore 2.5 mm, DIN, Material SBW	1 pz.
LLG09115613	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, with capillary side arms, NS 14.5, Bore 2.5 mm, DIN 12541-1, Material	1 pz.
LLG09115642	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 12.5, Bore 2.5 mm, DIN, Material	1 pz.
LLG06093232	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 14.5, Bore 4 mm, DIN, Material SBW	1 pz.
LLG06800209	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 18.8, Bore 4 mm, DIN, Material SBW	1 pz.
LLG09115614	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 14.5, Bore 4 mm, DIN, Material	1 pz.
LLG09115644	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 18.8, Bore 4 mm, DIN 12541-1, Material	1 pz.
LLG07510838	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 21.5*, Bore 6 mm, DIN, Material DURAN®	1 pz.
LLG09115615	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 18.8, Bore 6 mm, DIN, Material	1 pz.
LLG09115617	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 21.5*, Bore 6 mm, DIN 12541-1, Material	1 pz.

Code	Description	Packaging
LLG07657852	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 24, Bore 8 mm, DIN, Material DURAN®	1 pz.
LLG09115618	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 24*, Bore 8 mm, DIN, Material	1 pz.
LLG06266986	Stopcocks, with glass plug, borosilicate glass 3.3, Spare key, solid glass, NS 29.2, Bore 10 mm, DIN, Material DURAN®	1 pz.
LLG09115619	Stopcocks, with glass plug, borosilicate glass 3.3, With solid glass key, NS 29.2*, Bore 10 mm, DIN 12541-1, Material	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

