


**Stopcocks, with PTFE plug, borosilicate glass 3.3****Lenz****Stopcocks, with PTFE plug, borosilicate glass 3.3**

Acc. to DIN 12541-1 and in addition to DIN standard, with PTFE plug and screwthread retaining device, one straight bore, material: borosilicate glass 3.3, leakage rate acc. to DIN 12540-1.



Code	Description	Packaging
LLG07072146	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 12.5, Bore 1.5 mm, DIN, Material PTFE	1 pz.
LLG09115631	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 12.5, Bore 1.5 mm, DIN DIN 12541-1, Material	1 pz.
LLG09115632	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 12.5*, Bore 1.5 mm, DIN DIN 12541-1, Material	1 pz.
LLG06070033	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 12.5, Bore 2.5 mm, DIN, Material	1 pz.
LLG06306550	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 12.5, Bore 2.5 mm, DIN, Material PTFE	1 pz.
LLG06801361	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 14.5, Bore 2.5 mm, DIN, Material PTFE	1 pz.
LLG09115633	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 14.5*, Bore 2.5 mm, DIN DIN 12541-1, Material	1 pz.
LLG09115648	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 14.5, Bore 2.5 mm, DIN DIN 12541-1, Material	1 pz.
LLG06306551	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 14.5, Bore 4 mm, DIN, Material PTFE	1 pz.
LLG07018955	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 18.8, Bore 4 mm, DIN, Material PTFE	1 pz.
LLG09115634	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 14.5, Bore 4 mm, DIN, Material	1 pz.
LLG09115635	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 18.8, Bore 4 mm, DIN DIN 12541-1, Material	1 pz.
LLG06222868	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 18.8, Bore 6 mm, DIN, Material	1 pz.
LLG07630235	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 18.8, Bore 6 mm, DIN, Material PTFE	1 pz.
LLG09115637	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 21.5, Bore 6 mm, DIN DIN 12541-1, Material	1 pz.
LLG07657006	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 24, Bore 8 mm, DIN, Material PTFE	1 pz.
LLG09115638	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 24, Bore 8 mm, DIN, Material	1 pz.
LLG06266985	Stopcocks, with PTFE plug, borosilicate glass 3.3, Spare key, NS 29.2, Bore 10 mm, DIN, Material PTFE	1 pz.



Code	Description	Packaging
LLG09115639	Stopcocks, with PTFE plug, borosilicate glass 3.3, With PTFE key, NS 29.2, Bore 10 mm, DIN DIN 12541-1, Material	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

