

**Sampling lance milky, stainless steel V4A, acc. to EN ISO 707****Bürkle****Sampling lance milky, stainless steel V4A, acc. to EN ISO 707**

Universal powder and granules sampler collector made of AISI 316 stainless steel. Developed especially for milk powder samples, conforms to standards of DIN and IDF (Int. Dairy Federation). The collector is inserted into bags or casks and the sample is removed using a semi-circular profile lance. Autoclavable.

Lock seal close-it and close-it food & pharma please order separately.



Code	Description	Packaging
LLG06051629	Sampling lance milky, stainless steel V4A, acc. to EN ISO 707, 150 ml, Length 48 cm, Imm. depth 385 mm, Ext. diam. 28 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

