



Zone samplers, Multi/All-layer/Jumbo/Uno, steel V4A or aluminium

Bürkle

Zone samplers, Multi/All-layer/Jumbo/Uno, steel V4A or aluminium

For sampling bulk goods of all types. Available in various different lengths and materials:

Multi-sampler: With several closed chambers for multi-point sampling from several predetermined sampling depths. Available in stainless steel V4A/PTFE, stainless steel V4A/V4A or aluminium. Tube diameter: 25mm.

All-layer-sampler: with open inner tubes for cross-sectional sampling through all layers of the sample. Emptied via the open tube end of the sampler. Available in stainless steel V4A or aluminium. Tube diameter: 25mm.

Jumbo-sampler: generally as the All-layer-sampler with open tube end. With a tube diameter of 50 mm, ideal for coarse grains and larger quantities. For multi point sampling the chambers can be separated by PTFE separators.

Uno-sampler: With a closed chamber, for target samples from a pre-determined sample depth.

Easy emptying with the emptying device.

Lock seal close-it and close-it food & pharma please order separately.



Code	Description	Packaging
LLG09303886	Uno-Sampler, Chamber Capacity 17 ml, Length 55 cm, Imm. depth 430 mm, Number openings 1, in Stainless steel (V4A)/PTFE	1 pz.
LLG09303887	Uno-Sampler, Chamber Capacity 17 ml, Length 85 cm, Imm. depth 710 mm, Number openings 1, in Stainless steel (V4A)/PTFE	1 pz.
LLG09303888	Uno-Sampler, Chamber Capacity 17 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 1, in Stainless steel (V4A)/PTFE	1 pz.
LLG07631201	Multi-Sampler, Chamber Capacity 14 ml, Length 55 cm, Imm. depth 430 mm, Number openings 3, in V4A/V4A	1 pz.
LLG09303860	Multi-Sampler, Chamber Capacity 14 ml, Length 55 cm, Imm. depth 430 mm, Number openings 3, in Stainless steel (V4A)/PTFE	1 pz.
LLG09303866	Multi-Sampler, Chamber Capacity 17 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 5, in V4A/V4A	1 pz.
LLG09303870	Multi-Sampler, Chamber Capacity 14 ml, Length 55 cm, Imm. depth 430 mm, Number openings 3, in Aluminium	1 pz.
LLG07653124	Multi-Sampler, Chamber Capacity 17 ml, Length 85 cm, Imm. depth 710 mm, Number openings 3, in V4A/V4A	1 pz.
LLG09303862	Multi-Sampler, Chamber Capacity 17 ml, Length 85 cm, Imm. depth 710 mm, Number openings 3, in Stainless steel (V4A)/PTFE	1 pz.
LLG09303864	Multi-Sampler, Chamber Capacity 17 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 5, in Stainless steel (V4A)/PTFE	1 pz.
LLG09303872	Multi-Sampler, Chamber Capacity 17 ml, Length 85 cm, Imm. depth 710 mm, Number openings 3, in Aluminium	1 pz.
LLG09303874	Multi-Sampler, Chamber Capacity 17 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 5, in Aluminium	1 pz.
LLG06232190	Cross-section sampler, Chamber Capacity 70 ml, Length 55 cm, Imm. depth 430 mm, Number openings 3, in V4A/V4A	1 pz.
LLG09303875	Cross-section sampler, Chamber Capacity 70 ml, Length 55 cm, Imm. depth 430 mm, Number openings 3, in Aluminium	1 pz.
LLG06234598	Cross-section sampler, Chamber Capacity 130 ml, Length 85 cm, Imm. depth 710 mm, Number openings 3, in V4A/V4A	1 pz.

Code	Description	Packaging
LLG09303876	Cross-section sampler, Chamber Capacity 130 ml, Length 85 cm, Imm. depth 710 mm, Number openings 3, in Aluminium	1 pz.
LLG09303865	Cross-section sampler, Chamber Capacity 260 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 5, in V4A/V4A	1 pz.
LLG09303877	Cross-section sampler, Chamber Capacity 260 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 5, in Aluminium	1 pz.
LLG06240956	Jumbo sampler, Chamber Capacity 390 ml, Length 85 cm, Imm. depth 710 mm, Number openings 3, in Aluminium/PTFE	1 pz.
LLG07625884	Jumbo sampler, Chamber Capacity 650 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 5, in Aluminium/PTFE	1 pz.
LLG09303878	Jumbo sampler, Chamber Capacity 880 ml, Length 85 cm, Imm. depth 710 mm, Number openings 3, in Aluminium	1 pz.
LLG06236850	Jumbo sampler, Chamber Capacity 910 ml, Length 250 cm, Imm. depth 2355 mm, Number openings 7, in Aluminium/PTFE	1 pz.
LLG09303879	Jumbo sampler, Chamber Capacity 1700 ml, Length 150 cm, Imm. depth 1355 mm, Number openings 5, in Aluminium	1 pz.
LLG09303883	Jumbo sampler, Chamber Capacity 2900 ml, Length 250 cm, Imm. depth 2355 mm, Number openings 7, in Aluminium	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



