

**Sampler ProfiSampler, aluminium****Bürkle****Sampler ProfiSampler, aluminium**

Vacuum-operated sample collector for sampling solvents and flammable liquids (Explosion Categories IIA, IIB and IIC) out of drums, tanks, road tankers, railway wagons, ships, lakes, ground water.

- Filling level control, automatic return flow and overfill protection by built-in venting valve.
- Suction volume: 200 ml, suction pressure 0.5 bar and suction depth down to 4m of water.
- With bottle thread GL 45 (suitable for Schott glass bottles 100 to 1000ml)
- Can be supplied with optimized, straight suction tube for simpler sampling.
- Transport case for tidy transport and disposal
- Always check the chemical resistance of parts coming into contact with the sample.
- When taking samples of flammable liquids, the device must be earthed.

Wear protective clothing!

Scope of delivery: aluminium vacuum pump, with feet (aluminium), venting valve and GL45 bottle thread, 1 x 500 ml glass bottle, 5 x 220 mm long, conductive, black PA suction tubes, 1 tubing cutter, 1 brass hose weight, 20 seal-it bottle seals, earthing cable.



Code	Description	Packaging
LLG06225622	Sampler ProfiSampler, aluminium, Type ProfiSampler Aluminium	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

