



Accessories for Rotary Evaporators Hei-VAP Series

Accessories for Rotary Evaporators Hei-VAP Series



Code	Description	Packaging
LLG06256590	Accessories for Rotary Evaporators Hei-VAP Series, AquaStop - security mechanism for water disconnection	1 pz.
LLG07651893	Accessories for Rotary Evaporators Hei-VAP Series, Boiling temperature sensor	1 pz.
LLG06236805	Accessories for Rotary Evaporators Hei-VAP Series, Guard hood for evaporating flasks	1 pz.
LLG09812476	Accessories for Rotary Evaporators Hei-VAP Series, Guard shield for attaching to bath	1 pz.
LLG07055667	Accessories for Rotary Evaporators Hei-VAP Series, Spare PTFE vacuum seal	1 pz.
LLG09812480	Accessories for Rotary Evaporators Hei-VAP Series, Switchbox, connects up to 3 rotary evaporators to one pump and control vacuum individually	1 pz.
LLG09812479	Accessories for Rotary Evaporators Hei-VAP Series, Temperature sensor for controller function AUTOaccurate	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

