

**Eyewash Emergency Station, Wall-Mounting****B-SAFETY****Eyewash Emergency Station, Wall-Mounting**

A dust-proof wall mounted station containing one bottle of 200 ml sterile phosphate buffer solution and one bottle of 500 ml sterile sodium chloride solution complete with instructions and a mirror or eye wash emergency station DUO with one bottle of 500 ml sterile phosphate buffer solution and one bottle of 1000 ml sterile sodium chloride solution. Especially suitable for workplaces where on the spot rinsing is required, smaller workplaces. Bottles according to DIN EN 15154-4.

Dimensions eye wash station (W x D x H):

Dimensions eye wash station DUO (W x D x H): 227 x 110 x 265mm  
250 x 115 x 355mm

Durability: 3 years

Flush time: about 7 minutes



Code	Description	Packaging
LLG09733787	Eyewash Emergency Station, Wall-Mounting, Eyewash emergency station, wall-mounting	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

