

**Pigment Lance****Bürkle****Pigment Lance**

For pigments and pastes and agglomerating bulk goods.

- With half-open flute 50% for heavily stiching products
- With one-quarter-open flute 25% for lightly stiching products
- Stainless steel V4A (1.4404)



Code	Description	Packaging
LLG09303710	Pigment Lance, 145 ml, Flute opening 25%, Length 30 cm, Diam. 30 mm	1 pz.
LLG09303713	Pigment Lance, 179 ml, Flute opening 50%, Length 60 cm, Diam. 30 mm	1 pz.
LLG09303715	Pigment Lance, 273 ml, Flute opening 50%, Length 90 cm, Diam. 30 mm	1 pz.
LLG09303712	Pigment Lance, 312 ml, Flute opening 25%, Length 60 cm, Diam. 30 mm	1 pz.
LLG09303714	Pigment Lance, 480 ml, Flute opening 25%, Length 90 cm, Diam. 30 mm	1 pz.
LLG09303711	Pigment Lance, 85 ml, Flute opening 50%, Length 30 cm, Diam. 30 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

