

**3-channel-short period timer****Dostmann electronic****3-channel-short period timer**

Indispensable for laboratories. Triple timer with large display.

- Indicates three separate times simultaneously
- Measuring range: 1 sec to 20 hours
- Stopwatch function
- Integrated quartz watch
- Supplied with magnet and clip
- Operating temperature: 0 to 50 °C
- Battery 1.5 V AAA inclusive

Dimensions (W x H x D): 71 x 88 x 24mm

Weight: 115g



| Code        | Description   | Packaging |
|-------------|---|-----------|
| LLG06233690 | 3-channel timer, Measuring capacity 20:00:00 hr: min: sec, Acoustic alarm Yes | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

