



Control units bioMIXcontrol / bioMIXcontrol S for stirring drives bioMIXdrive

The powerful control units are characterized by:

- Fully encapsulated and sealed stainless steel housing, water, dust and germ-proof
- SoftStart for reliable catching / centering and safe acceleration of the stirring rod with different volumes and viscosities
- Jerk-free stirring even at low speeds
- QuickSet for quick setting of start- and maximum speed
- QuickStop-function against after-run and glass breakage, storage of last working status when turning off
- Digital display for speed setting
- Extremely robust construction and performance
- Easy to clean and maximum resistance
- 3 years warranty
- Made in Germany

Power setting (10-steps) for high power with regard to large and viscous stirring amounts and reduced power for the warming-free continuous operation of stirring drives bioMIXdrive 1, 2, 3 or 4. The model bioMIXcontrol S is also stackable.

Specifications

Speed: 5 ... 250 rpm

Stirring power: 40 W

Power setting: 10 ... 100 % (10 steps)

Weight: approx. 1.3 kg



Code	Description	Packaging
LLG09645982	Control unit bioMixcontrol 4MS for bioMIXdrive 1 only	1 pz.
LLG09645951	Control unit bioMIXcontrol 5-250 rpm, 10-stepped	1 pz.
LLG09645981	Control unit bioMixcontrol S stapable	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

