

**Magnetic stirrers for cell cultures with external control, bioMIXdrive 1/2/3/4**

Due to the external control, these magnetic stirrers are suitable for use in ovens, incubators, safety cabinets or for underwater operation, depending on the device. Wide speed ranges make them ideal for viscous media. Low space requirement due to the flat and compact design with maximum stirring surface.

- Maintenance-free 2mag motor drive
- Fully sealed housing
- Extremely robust construction and design
- Easy to clean and maximum resistance
- 3 years warranty
- Made in Germany

**Please always order control units separately.**

**Submersible magnetic stirrers for careful, protective and warming-free mixing of cell cultures and culture broths.**

- Suitable for propeller and spinning ball culture flasks
- For long-term use, e.g., in CO<sub>2</sub> incubators
- Individual measurements on request

**Specifications**

**Stirring quantity (H<sub>2</sub>O): 5 ... 5000 ml**

**Permissible ambient temperature: -10 up to +50 °C at 100 % H.R. max.**

**+50 °C submerged in water**

**Housing: Stainless steel**

**Protection class: IP 68**



Code	Description	Packaging
LLG09645979	Cell culture stirrer bioMixdrive 3 5-250 rpm	1 pz.
LLG09645950	Cellculture-stirrer bioMixdrive 1 5-250 rpm	1 pz.
LLG09645978	Cell culture stirrer bioMixdrive 2 5-250 rpm	1 pz.
LLG09645980	Cell culture stirrer bioMixdrive 4 5-250 rpm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

