

**LLG-Funnel holder, PP****LLG Labware****LLG-Funnel holder, PP**

Manufactured from best quality polypropylene. Suitable to hold all types of laboratory funnels with diameters up to 150mm. Two tapered wedges are provided to ensure firm grip on rods having diameters of 10-13mm.



Code	Description	Packaging
LLG09251558	Funnel holder, PP, To hold funnels 1 piece	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

