



Inlet filter for solvents

S.C.A.T.

Inlet filter for solvents

HPLC solvent filters for 3.2mm OD tubing.



Code	Description	Packaging
LLG04005891	Inlet filter for solvents, HPLC solvent filter, PFA/PTFE, for 1/8" (3.2mm o.d.)	5 pz.
LLG04005890	Inlet filter for solvents, HPLC solvent filter, PP, for 1/8" dia. (3.2mm o.d.)	5 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

