




LLG Septa for Screw Caps ND24

LLG Labware

LLG Septa for Screw Caps ND24



Code	Description	Packaging
LLG04008291	LLG Septa for Screw Caps ND24, Cap size ND24, Ø 22 mm, Material Silicone/PTFE, Colour White/beige, Shore hardness 45 A, Thickness 3.2 mm	100 pz
LLG09003542	LLG Septa for Screw Caps ND24, Cap size ND24, Ø 22 mm, Material Butyl/PTFE, Colour Red/grey, Shore hardness 55 A, Thickness 2.5 mm	1000 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

