




**LLG Septa for Screw Caps ND24**

**LLG Labware**

**LLG Septa for Screw Caps ND24**



| Code        | Description                                                                                                                              | Packaging |
|-------------|------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG04008291 | LLG Septa for Screw Caps ND24, Cap size ND24, Ø 22 mm, Material Silicone/PTFE, Colour White/beige, Shore hardness 45 A, Thickness 3.2 mm | 100 pz    |
| LLG09003542 | LLG Septa for Screw Caps ND24, Cap size ND24, Ø 22 mm, Material Butyl/PTFE, Colour Red/grey, Shore hardness 55 A, Thickness 2.5 mm       | 1000 pz   |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

