

**LLG-PP Screw Seals ND24, UltraBond****LLG Labware****LLG-PP Screw Seals ND24, UltraBond**

Cap and liner form an inseparable unit, so that the liner cannot be pushed into the vial even with a blunt needle.



Code	Description	Packaging
LLG04008292	LLG-PP Screw Seals ND24, UltraBond, Cap size ND24, Material PP, Colour White, Without hole, Material septa Silicone/PTFE, Colour septa White/beige, Thickness 3.20 mm, Shore hardness 45 A, Septa type UltraBond	100 pz
LLG04008293	LLG-PP Screw Seals ND24, UltraBond, Cap size ND24, Material PP, Colour White, With hole, Int. bore Ø 15.0 mm, Material septa Silicone/PTFE, Colour septa White/beige, Thickness 3.20 mm, Shore hardness 45 A, Septa type UltraBond	100 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

