



Rotors and Receptacles SM 100 / 200 / 300

Retsch

Rotors and Receptacles SM 100 / 200 / 300



Code	Description	Packaging
LLG09739074	Rotors and Receptacles SM 100 / 200 / 300, 6-disc rotor stainless steel, with reversible cutting tips	1 pz.
LLG09739607	Rotors and Receptacles SM 100 / 200 / 300, Collecting receptacle, 30L plastic, incl. adapter and filter hose	1 pz.
LLG09739608	Rotors and Receptacles SM 100 / 200 / 300, Collecting unit 250/500ml incl. 2 sample bottles	1 pz.
LLG07900427	Rotors and Receptacles SM 100 / 200 / 300, Parallel section rotor stainless steel	1 pz.
LLG06267310	Rotors and Receptacles SM 100 / 200 / 300, V rotor stainless steel for SM 300	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

