

**Accessories for Ultrasonic Processor UP200H****Hielscher Ultrasonics****Accessories for Ultrasonic Processor UP200H**

Sonotrode S3, titanium, approx. length 100mm

Sonotrode S7, titanium, approx. length 100mm

Sonotrode S14, titanium, approx. length 100mm

Sonotrode S40, titanium, approx. length 100mm

Stand ST1-16, rod diameter 16mm, stainless steel, base (L x W x H) 300 x 150 x 600mm

Sound Protection Box SB3-16, with adjustable table and 16mm stand pole

Timer T1, from 00:00 - 99:59 (min:sec)



Code	Description	Packaging
LLG04005992	Accessories for Ultrasonic Processor UP200H, Type S40, Tip dia. 40mm, For volume 100 - 2000 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

