

**Mini Plate Spinner MPS 1000****Corning BV****Mini Plate Spinner MPS 1000**

Centrifuge designed specifically for quick and easy centrifugation of samples in PCR plates or microplates.

- Capacity for two standard plates
- Accepts skirted, non-skirted and all standard PCR plates
- Use before and after thermal cycling to increase PCR yield
- Top loading design
- Exceptionally small footprint

Specifications

Max. RCF: 500xg

Max. speed: 2500rpm

Max. capacity: 2 PCR Plates

Power consumption: 75W

Dimensions (W x D x H): 160 x 160 x 180 mm

Weight: 1.5 kg

Power supply: 230V 50/60Hz



Code	Description	Packaging
LLG09945790	Mini Plate Spinner MPS 1000, Type MPS 1000	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

