

**Erlenmeyer flasks DURAN Super Duty, wide neck****DWK Life Sciences****Erlenmeyer flasks DURAN® Super Duty, wide neck**

With approximate graduations. Reinforced, sturdy rim and increased mechanical strength for heavy duty. Glass Type I / neutral glass as per USP, EP and JP. With retrace code (Batch Identification). Certificate available via the Internet. Autoclavable.



Code	Description	Packaging
LLG09141280	Erlenmeyer flasks DURAN® Super Duty, wide neck, 100 ml, Neck diam. 34 mm, Bottom diam. 64 mm, Height 105 mm	1 pz
LLG09141281	Erlenmeyer flasks DURAN® Super Duty, wide neck, 250 ml, Neck diam. 50 mm, Bottom diam. 85 mm, Height 140 mm	1 pz
LLG09141282	Erlenmeyer flasks DURAN® Super Duty, wide neck, 500 ml, Neck diam. 50 mm, Bottom diam. 105 mm, Height 175 mm	1 pz
LLG09141283	Erlenmeyer flasks DURAN® Super Duty, wide neck, 1000 ml, Neck diam. 50 mm, Bottom diam. 131 mm, Height 220 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

