

**Temperature / Humidity data logger LOG 110****Dostmann electronic****Temperature / Humidity data logger LOG 110**

Large display for actual value, battery status display, easy to use software. Perfect designed instrument for transportation or storage after HACCP.

- Quality surveillance in laboratories, at the production or in the storage
- Ambient conditions in manufacturing facilities
- Ambient conditions in power plants
- Monitoring of fridges for the pharma industry
- Memory for 60000 data sets

Please order Software with USB-Cabel separately (9.726 413).

Specifications**Temperature**

Measuring range int.:

Measuring range ext.: -30 ... 70 °C

-50 ... 125 °C

Accuracy: 0.1 °C

Precision: Int.: ± 0.5 °C (-10 ... 40 °C), otherwise ± 0.7 °C

Ext.: ± 0.5 °C (0 ... 40 °C), otherwise ± 0.7 °C

Humidity

Measuring range: 0 ... 99 % r.H.

Resolution: 0.1% r.H.

Accuracy: ± 3 % r.H.

Housing: Plastic (ABS)

Power supply: 1 x CR2032 3 V

Dimensions: 88 x 55 x 20 mm

Weight: 95 g (incl. packaging)

Alarm function: Alarm for too low or too high temperature and/or humidity values

Calibration reminder



Code	Description	Packaging
LLG09726412	Temperature / Humidity data logger LOG 110, Type LOG 110, Range start temperature -30 °C	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

