

**Tube Rack pop-up, foldable, PP**

Heathrow Scientific

**Tube Rack pop-up, foldable, PP**

Robust polypropylene construction. Simple, easy to open and close. Flatten for compact storage. Pops up for instant use. Expand and collapse when racks are needed. Printed grid aids rapid retrieval. Stack racks in open or closed position. Saves space when racks are not in use. Not recommended for steam sterilisation. Dimensions (L x W x H): 255mm x 137mm x 72mm when expanded, 21mm high when collapsed.



Code	Description	Packaging
LLG09193998	Tube Rack pop-up, foldable, PP, For 21 x 15ml and 12 x 50ml tubes, Colour Blue	2 pz.
LLG06254578	Tube Rack pop-up, foldable, PP, For 18 x 50 ml tubes, Colour Green	2 pz.
LLG09193997	Tube Rack pop-up, foldable, PP, For 21 x 15ml and 12 x 50ml tubes, Colour Pink	2 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

