



Round bottom flasks, borosilicate glass 3.3

Round bottom flasks, borosilicate glass 3.3

Resistant to heat and almost all chemicals. According to DIN EN ISO 4797.

Reinforced rims are used to prevent breaking during washing and transportation.



Code	Description	Packaging
LLG04008385	Round bottom flasks, borosilicate glass 3.3, 500 ml, Flask diam. 105 mm, Neck diam. approx. 30 mm, Height 175 mm, Socket 24/29 NS	1 pz.
LLG04008381	Round bottom flasks, borosilicate glass 3.3, 100 ml, Flask diam. 64 mm, Neck diam. approx. 22 mm, Height 115 mm, Socket 19/26 NS	1 pz.
LLG04008382	Round bottom flasks, borosilicate glass 3.3, 100 ml, Flask diam. 64 mm, Neck diam. approx. 30 mm, Height 115 mm, Socket 24/29 NS	1 pz.
LLG06236346	Round bottom flasks, borosilicate glass 3.3, 100 ml, Flask diam. 64 mm, Neck diam. approx. 22 mm, Height 115 mm, Socket 14/23 NS	1 pz.
LLG06236774	Round bottom flasks, borosilicate glass 3.3, 100 ml, Flask diam. 64 mm, Neck diam. approx. 34 mm, Height 115 mm, Socket 29/32 NS	1 pz.
LLG06237702	Round bottom flasks, borosilicate glass 3.3, 1000 ml, Flask diam. 131 mm, Neck diam. approx. 42 mm, Height 210 mm, Socket 24/29 NS	1 pz.
LLG06242240	Round bottom flasks, borosilicate glass 3.3, 1000 ml, Flask diam. 131 mm, Neck diam. approx. 42 mm, Height 210 mm, Socket 29/32 NS	1 pz.
LLG06238472	Round bottom flasks, borosilicate glass 3.3, 2000 ml, Flask diam. 165 mm, Neck diam. approx. 34 mm, Height 260 mm, Socket 29/32 NS	1 pz.
LLG04008383	Round bottom flasks, borosilicate glass 3.3, 250 ml, Flask diam. 85 mm, Neck diam. approx. 22 mm, Height 145 mm, Socket 19/26 NS	1 pz.
LLG04008384	Round bottom flasks, borosilicate glass 3.3, 250 ml, Flask diam. 85 mm, Neck diam. approx. 30 mm, Height 145 mm, Socket 24/29 NS	1 pz.
LLG06236303	Round bottom flasks, borosilicate glass 3.3, 250 ml, Flask diam. 85 mm, Neck diam. approx. 34 mm, Height 145 mm, Socket 29/32 NS	1 pz.
LLG04008379	Round bottom flasks, borosilicate glass 3.3, 50 ml, Flask diam. 51 mm, Neck diam. approx. 22 mm, Height 105 mm, Socket 19/26 NS	1 pz.
LLG04008380	Round bottom flasks, borosilicate glass 3.3, 50 ml, Flask diam. 51 mm, Neck diam. approx. 30 mm, Height 105 mm, Socket 24/29 NS	1 pz.
LLG06236508	Round bottom flasks, borosilicate glass 3.3, 50 ml, Flask diam. 51 mm, Neck diam. approx. 22 mm, Height 105 mm, Socket 14/23 NS	1 pz.

Code	Description	Packaging
LLG06238345	Round bottom flasks, borosilicate glass 3.3, 50 ml, Flask diam. 51 mm, Neck diam. approx. 34 mm, Height 105 mm, Socket 29/32 NS	1 pz.
LLG06240426	Round bottom flasks, borosilicate glass 3.3, 500 ml, Flask diam. 105 mm, Neck diam. approx. 34 mm, Height 175 mm, Socket 29/32 NS	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

