



Disposable Erlenmeyer Flasks Nalgene with vented closure, Type 4115, 4116, PETG, sterile

Thermo Scientific

Disposable Erlenmeyer Flasks Nalgene with vented closure, Type 4115, 4116, PETG, sterile

With blue vented closure, the chance for cross contamination is reduced. Ideal for shaker and suspension cell culture, media preparation or storage. Made of light, crystal clear PETG. Molded-in graduations. Leak-proof HDPE closure has a hydrophobic 0.2µm PTFE membrane that allows sterile air exchange. Flasks offer a 5-year shelf life, a 10⁻⁶ SAL, are non-pyrogenic and non-cytotoxic. Available with flat or baffled bottom. Individually packaged for easy storage and handling.



Code	Description	Packaging
LLG09141305	Disposable Erlenmeyer Flasks Type 4115, Flat Bottom, 1000 ml, Cover type 45-430	6 pz.
LLG09141311	Disposable Erlenmeyer Flasks Type 4116, Baffled Bottom, 1000 ml, Cover type 45-430	6 pz.
LLG06234530	Disposable Erlenmeyer Flasks Type 4115, Flat Bottom, 125 ml, Cover type 38-430	24 pz.
LLG09141308	Disposable Erlenmeyer Flasks Type 4116, Baffled Bottom, 125 ml, Cover type 38-430	24 pz.
LLG09141306	Disposable Erlenmeyer Flasks Type 4115, Flat Bottom, 2000 ml, Cover type 45-430	4 pz.
LLG09141312	Disposable Erlenmeyer Flasks Type 4116, Baffled Bottom, 2000 ml, Cover type 45-430	4 pz.
LLG06231087	Disposable Erlenmeyer Flasks Type 4115, Flat Bottom, 250 ml, Cover type 38-430	12 pz.
LLG09141309	Disposable Erlenmeyer Flasks Type 4116, Baffled Bottom, 250 ml, Cover type 38-430	12 pz.
LLG09141307	Disposable Erlenmeyer Flasks Type 4115, Flat Bottom, 2800 ml, Cover type 70	4 pz.
LLG09141313	Disposable Erlenmeyer Flasks Type 4116, Baffled Bottom, 2800 ml, Cover type 70	4 pz.
LLG09141304	Disposable Erlenmeyer Flasks Type 4115, Flat Bottom, 500 ml, Cover type 45-430	12 pz.
LLG09141310	Disposable Erlenmeyer Flasks Type 4116, Baffled Bottom, 500 ml, Cover type 45-430	12 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

