



Stopcocks LaboPlast

Bürkle

Stopcocks LaboPlast®**StopCock Siphon**

- With rotatable spout, preventing drips
- With tapered spout for bottles with narrow openings from diameter 6mm to 8mm
- Fits screw threads with G $\frac{3}{4}$ "
- Pressure proof to 2 bar at 20°C (without discharge bracket and grommet)
- Parts in contact with media made of PP, PE and PTFE
- Food safe

StopCock Compact

- Compact design, the short discharge does not protrude out over the lower edge of the container
- The specially optimised design of the discharge forms an even stream of liquid with no splashing or spattering
- The Compact stop cock fits all containers with G $\frac{3}{4}$ " connecting threads
- PP, PE and PTFE are materials that offer outstanding chemical resistance to acids and bases
- Food safe

StopCock PA

- Made of polyamide, electroconductive, can be grounded separately
- Ideally suited for solvents
- Inexpensive alternative to stainless steel stopcocks
- PTFE seal
- Including adapter for $\frac{3}{4}$ " inner thread made of PA electroconductive
- Fits all containers with $\frac{3}{4}$ " outer thread and with included adapter, to all $\frac{3}{4}$ " inner threads



Code	Description	Packaging
LLG09139997	Stopcocks LaboPlast®, Type StopCock Siphon	1 pz
LLG09140000	Stopcocks LaboPlast®, Type StopCock PA	1 pz
LLG09140003	Stopcocks LaboPlast®, Type StopCock Compact	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

