



## Stopcocks LaboPlast

Bürkle

**Stopcocks LaboPlast®****StopCock Siphon**

- With rotatable spout, preventing drips
- With tapered spout for bottles with narrow openings from diameter 6mm to 8mm
- Fits screw threads with G  $\frac{3}{4}$ "
- Pressure proof to 2 bar at 20°C (without discharge bracket and grommet)
- Parts in contact with media made of PP, PE and PTFE
- Food safe

**StopCock Compact**

- Compact design, the short discharge does not protrude out over the lower edge of the container
- The specially optimised design of the discharge forms an even stream of liquid with no splashing or spattering
- The Compact stop cock fits all containers with G  $\frac{3}{4}$ " connecting threads
- PP, PE and PTFE are materials that offer outstanding chemical resistance to acids and bases
- Food safe

**StopCock PA**

- Made of polyamide, electroconductive, can be grounded separately
- Ideally suited for solvents
- Inexpensive alternative to stainless steel stopcocks
- PTFE seal
- Including adapter for  $\frac{3}{4}$ " inner thread made of PA electroconductive
- Fits all containers with  $\frac{3}{4}$ " outer thread and with included adapter, to all  $\frac{3}{4}$ " inner threads



Code	Description	Packaging
LLG09140003	Stopcocks LaboPlast®, Type StopCock Compact	1 pz.
LLG09140000	Stopcocks LaboPlast®, Type StopCock PA	1 pz.
LLG09139997	Stopcocks LaboPlast®, Type StopCock Siphon	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

