

**LLG-Griffin beakers, PP****LLG Labware****LLG-Griffin beakers, PP**

Transparent. Low form with practical spout and easy-to-read blue graduation.

- Autoclavable at 121 °C
- Good chemical resistance
- Stackable
- According to ISO 7056



| Code        | Description                                                 | Packaging |
|-------------|-------------------------------------------------------------|-----------|
| LLG09013675 | Griffin beakers, PP, Nominal capacity 1000 ml, Grad. 50 ml  | 2 pz      |
| LLG09013676 | Griffin beakers, PP, Nominal capacity 2000 ml, Grad. 100 ml | 1 pz      |
| LLG09013670 | Griffin beakers, PP, Nominal capacity 25 ml, Grad. 5 ml     | 2 pz      |
| LLG09013671 | Griffin beakers, PP, Nominal capacity 50 ml, Grad. 10 ml    | 2 pz      |
| LLG09013672 | Griffin beakers, PP, Nominal capacity 100 ml, Grad. 10 ml   | 2 pz      |
| LLG09013674 | Griffin beakers, PP, Nominal capacity 500 ml, Grad. 50 ml   | 2 pz      |
| LLG09013673 | Griffin beakers, PP, Nominal capacity 250 ml, Grad. 25 ml   | 2 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

