

**LLG-Griffin beakers, PP****LLG Labware****LLG-Griffin beakers, PP**

Transparent. Low form with practical spout and easy-to-read blue graduation.

- Autoclavable at 121 °C
- Good chemical resistance
- Stackable
- According to ISO 7056



Code	Description	Packaging
LLG09013675	Griffin beakers, PP, Nominal capacity 1000 ml, Grad. 50 ml	2 pz
LLG09013676	Griffin beakers, PP, Nominal capacity 2000 ml, Grad. 100 ml	1 pz
LLG09013670	Griffin beakers, PP, Nominal capacity 25 ml, Grad. 5 ml	2 pz
LLG09013671	Griffin beakers, PP, Nominal capacity 50 ml, Grad. 10 ml	2 pz
LLG09013672	Griffin beakers, PP, Nominal capacity 100 ml, Grad. 10 ml	2 pz
LLG09013674	Griffin beakers, PP, Nominal capacity 500 ml, Grad. 50 ml	2 pz
LLG09013673	Griffin beakers, PP, Nominal capacity 250 ml, Grad. 25 ml	2 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

