




## Magnetic Stirring Beakerliner, PTFE

### Magnetic Stirring Beakerliner, PTFE

PTFE-encapsulated magnetic core (Alnico 5) axially mounted in a guide cage made of PTFE, universal chemical resistance. No shear action on the bottom of the beaker, smooth running in glass beakers also on an uneven bottom. The cage acts like a baffel and thus provides optimum mixing results.



| Code        | Description                                                                                                                | Packaging |
|-------------|----------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG09197587 | Magnetic Stirring Beakerliner, PTFE, For beaker 2000 ml, Ext. diam. 103 mm, Height 32 mm, Length 80 mm, Diam. 10 mm        | 1 pz.     |
| LLG09197588 | Magnetic Stirring Beakerliner, PTFE, For beaker 3000, 5000 ml, Ext. diam. 125 mm, Height 48 mm, Length 106 mm, Diam. 25 mm | 1 pz.     |
| LLG09197585 | Magnetic Stirring Beakerliner, PTFE, For beaker 400 ml, Ext. diam. 67 mm, Height 21 mm, Length 50 mm, Diam. 8 mm           | 1 pz.     |
| LLG09197586 | Magnetic Stirring Beakerliner, PTFE, For beaker 600 ml, Ext. diam. 74 mm, Height 28 mm, Length 60 mm, Diam. 9 mm           | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

