

**Precision balance Type PCB****Kern & Sohn****Precision balance Type PCB**

Compact scale for your laboratory. Due to the Recipe and piece counting function, as well as the grand precision this scale is perfect for the daily use in the Laboratory.

- Protective working cover over keyboard and housing is enclosed o standard
 - Modern, backlit LCD display, digit height 15 mm
 - Optional battery operation, 9V Block, not included
 - AUTO-OFF funktion to preserve the battery, can be switched off
 - Stainless steel weighing plate, the models PCB 100-3, PCB 250-3 and PCB 350-3 has a plastic weighing plate
- Overall Dimensions (W x D x H): 163 x 245 x 85mm
Weight: approx. 1kg



Code	Description	Packaging
LLG07619757	Precision balance Type PCB, Type PCB 100-3, 100 g	1 pz.
LLG06228513	Precision balance Type PCB, Type PCB 1000-2, 1000 g	1 pz.
LLG06235098	Precision balance Type PCB, Type PCB 1000-1, 1000 g	1 pz.
LLG07625237	Precision balance Type PCB, Type PCB 10000-1, 10000 g	1 pz.
LLG07635868	Precision balance Type PCB, Type PCB 200-2, 200 g	1 pz.
LLG06229399	Precision balance Type PCB, Type PCB 2000-1, 2000 g	1 pz.
LLG06235556	Precision balance Type PCB, Type PCB 2500-2, 2500 g	1 pz.
LLG06241856	Precision balance Type PCB, Type PCB 350-3, 350 g	1 pz.
LLG06240408	Precision balance Type PCB, Type PCB 3500-2, 3500 g	1 pz.
LLG09901300	Precision balance Type PCB, Type PCB 6000-0, 6000 g	1 pz.
LLG06231685	Precision balance Type PCB, Type PCB 6000-1, 6000 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

