



Hanging Scales HDB-N

Hanging Scales HDB-N

Convenient for hanging loads.

- Ideal for quick checking in goods receiving and shipping, also suitable for cases
- Data hold function: For easy reading of the weight, the display can be "frozen" in different ways. Either automatically when the weighing value remains unchanged or manually by pressing the Hold key
- Hook (stainless steel), can be hinged
- Batteries included, 2 x 1.5 V AAA, operating time approx. 180 h, AUTO OFF function to preserve the battery (can be switched off)



Code	Description	Packaging
LLG06231173	Hanging Scales HDB-N, Type HDB 5K5N, 5 kg, Readability 5 g, Weight 200 g, dim. (W x D x H) 72 x 25 x 105 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

