


**Magnetic stirring bars ASTEROID, triangular, PTFE****43953****Magnetic stirring bars ASTEROID, triangular, PTFE**

More than 4x higher torque transmissible. Unbeatably enhanced mixing effect in comparison to conventional stirring bars. Strong, stable, long-lasting magnetic field, with no demagnetization effect. Concave design for minimal surface contact and reliable rotation, triangular shape for powerful shear force and enhanced mixing at lower rotational speeds. Significantly higher rotation speeds achievable, high quality, even and sealed PTFE encapsulation. Conform to FDA requirements, steam sterilisable at 121 °C. Multiple enhancements against all types of magnetic stirrers.



Code	Description	Packaging
LLG09646057	Magnetic stirring bars ASTEROID, triangular, PTFE, Diam. 12 mm, Length 25 mm	1 pz
LLG09646040	Magnetic stirring bars ASTEROID, triangular, PTFE, Diam. 32 mm, Length 40 mm	1 pz
LLG09646058	Magnetic stirring bars ASTEROID, triangular, PTFE, Diam. 39 mm, Length 70 mm	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

