

**Multi-position Magnetic stirrer with internal control MIX 4 MS****43953****Multi-position Magnetic stirrer with internal control MIX 4 MS**

Powerful 4-position stirring, 100 % maintenance-free and wear-free, 2mag magnetic induction drive concept, **with individual and independent speed setting for each stirring position** (Master-Slave-Technology). On-/Off function, synchronous speed mode for all positions, extremely wide speed range of 100 to 2000 min⁻¹, jerk-free stirring even at low speeds. Clear digital display for settings of stirrer speeds, SoftStart procedure for reliable stirrer bar coupling/centring and safe acceleration, automatic memory function to save the last settings, outstanding long working life of the hermetically sealed drive system and resistant stainless steel housing for easy and fast cleaning.

Specifications**Max. stirring volume (H₂O): 3 l****Motor rating input: 28 W****Stirring positions: 4****Speed range: 100 ... 2000 min⁻¹****Permissible ambient temperature: -10 ... 50 °C****Permissible humidity: 80 %****Display: digital****Dimensions (W x D x H): 245 x 425 x 38 mm****Weight: 7,8 kg****IP code: IP 65****Power supply: 115/230 V, 50/60 Hz**

Code	Description	Packaging
LLG09646010	Multiple Magnetic stirrer MIX 4 MS 4 independent stirr places	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

