



LLG-General purpose thermometer, red filling

LLG Labware

LLG-General purpose thermometer, red filling

Graduation resistant against acids and alkalies.



Code	Description	Packaging
LLG09236720	Thermometer, general purpose -35...	1 pz
LLG09236723	Thermometer, general purpose -20...	1 pz
LLG09236725	Thermometer, general purpose -10...	1 pz
LLG09236730	Thermometer, general purpose -10...	1 pz.
LLG09236731	LLG Thermometer, general purpose -20...	1 pz
LLG09236721	Thermometer, general purpose -10..	1 pz
LLG09236722	Thermometer, general purpose -10...	1 pz
LLG09236724	Thermometer, general purpose -10...	1 pz
LLG09236727	Thermometer, general purpose -10...	1 pz
LLG09236726	Thermometer, general purpose -10...	1 pz
LLG06206594	General purpose thermometer, red filling, Measuring range -20 ... 110 °C, Graduations 1 °C, Length 300 mm, Diam. 6 ±1 mm, Immersion depth Total , Capillary White enameled, round	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

