

**LLG-Cryogenic storage boxes, plastic coated, 136 x 136****LLG Labware****LLG-Cryogenic storage boxes, plastic coated, 136 x 136**

Freezer and cryogenic storage boxes (CryoBoxes). Dimensions: 136 x 136mm square. For the storage of samples. Made of cardboard, **plastic coated** without divider. With lid. Further colours available on request.



Code	Description	Packaging
LLG04654548	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Blue, Int. height 130 mm	
LLG06077875	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Blue, Int. height 50 mm	1 pz
LLG06260006	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Blue, Int. height 75 mm	1 pz
LLG06260010	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Blue, Int. height 100 mm	1 pz
LLG06269788	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Blue, Int. height 32 mm	1 pz.
LLG06081021	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Green, Int. height 50 mm	1 pz
LLG06269790	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Green, Int. height 32 mm	1 pz.
LLG04654549	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Red, Int. height 130 mm	
LLG06081022	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Red, Int. height 50 mm	1 pz
LLG06260007	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Red, Int. height 75 mm	1 pz
LLG06260011	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Red, Int. height 100 mm	1 pz
LLG06269789	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Red, Int. height 32 mm	1 pz.
LLG04653682	Cryogenic storage boxes, plastic coated, 136 x 136, Colour White, Int. height 32 mm	1 pezzo
LLG06080130	Cryogenic storage boxes, plastic coated, 136 x 136, Colour White, Int. height 50 mm	1 pz
LLG06083519	Cryogenic storage boxes, plastic coated, 136 x 136, Colour White, Int. height 130 mm	1 pz
LLG06237840	Cryogenic storage boxes, plastic coated, 136 x 136, Colour White, Int. height 100 mm	1 pz.
LLG06802726	Cryogenic storage boxes, plastic coated, 136 x 136, Colour White, Int. height 75 mm	1 pz
LLG04654547	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Yellow, Int. height 100 mm	1 pz.
LLG06077876	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Yellow, Int. height 50 mm	1 pz
LLG06260009	Cryogenic storage boxes, plastic coated, 136 x 136, Colour Yellow, Int. height 75 mm	1 pz





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

