



LLG-Cryogenic storage boxes, plastic coated, 133 x 133

LLG-Cryogenic storage boxes, plastic coated, 133 x 133

Dimensions: 133 x 133mm square. For the storage of samples. Made of cardboard, **plastic coated**. Without divider. With lid. Further sizes and colours available on request.



Code	Description	Packaging
LLG06700568	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Blue, Int. height 50 mm	1 pz
LLG09698776	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Blue, Int. height 75 mm	1 pz.
LLG09698786	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Blue, Int. height 100 mm	1 pz.
LLG04654505	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Green, Int. height 32 mm	1 pz.
LLG04654509	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Green, Int. height 130 mm	
LLG06084538	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Green, Int. height 50 mm	1 pz
LLG06803056	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Green, Int. height 75 mm	1 pz
LLG09698788	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Green, Int. height 100 mm	1 pz.
LLG06081654	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Red, Int. height 50 mm	1 pz
LLG06260030	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Red, Int. height 75 mm	1 pz
LLG09698787	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Red, Int. height 100 mm	1 pz.
LLG06082701	Cryogenic storage boxes, plastic coated, 133 x 133, Colour White, Int. height 50 mm	1 pz
LLG06260029	Cryogenic storage boxes, plastic coated, 133 x 133, Colour White, Int. height 32 mm	1 pz
LLG06260032	Cryogenic storage boxes, plastic coated, 133 x 133, Colour White, Int. height 100 mm	1 pz
LLG06260033	Cryogenic storage boxes, plastic coated, 133 x 133, Colour White, Int. height 130 mm	1 pz
LLG06801707	Cryogenic storage boxes, plastic coated, 133 x 133, Colour White, Int. height 75 mm	1 pz
LLG06260031	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Yellow, Int. height 75 mm	1 pz
LLG06260091	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Yellow, Int. height 50 mm	1 pz.
LLG09698789	Cryogenic storage boxes, plastic coated, 133 x 133, Colour Yellow, Int. height 100 mm	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

