



LLG-PE Caps ND20, transparent, for DIN crimp neck

LLG Labware

LLG-PE Caps ND20, transparent, for DIN crimp neck



Code	Description	Packaging
LLG07052184	LLG-PE Caps ND20, transparent, for DIN crimp neck, With hole, Int. bore Ø 4.3 mm	100 pz
LLG09003545	LLG-PE Caps ND20, transparent, for DIN crimp neck, With hole, Int. bore Ø 4.3 mm, Material septa Butyl/PTFE, Colour septa Red/grey, Thickness 1.30 mm, Shore hardness 55 A	1000 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

